INSTRUCTION FOR USE
사용설명서

1. Manufacturer Name
제조사

DR-VU (Co.,Ltd)
(주)디알뷰

Registered no. 5882

2. Manufacturer Address Manufacturing site
주소

40, Dongil-ro 96-gil, Jungnang-gu,Seoul, Republic of Korea
서울시 중랑구 동일로96길40, B1

3. Product name
품목명

Lead-free Shield,radiation, hand
손방사선방어기구

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6. Intended purpose
사용목적

This is a tool used to defend against radiation exposure to the hands

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INSTRUCTION FOR USE
사용설명서

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A. Preparations before use

a. Before use, check the condition of the product to see if there is any damage before wearing it.

b. The users should study how to use it, if not knowing how to use it, get help from other staff.

B. How to use it

Use the product's Velcro to properly wear it on the user.

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1. Be careful not to scratch the surface.
2. Don't fold and leave it pressed by other heavy things.
3. Do not use it for any purpose other than the purpose.
4. Inspect periodically and replace any cracks or tears.

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WHEN REUSE

The device is supplied non-sterile and be reused after simple surface cleaning with water paper or alcohol paper if any stain.

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1. After using it, with wet tissue or soft cloth, wipe it off, and dry it sufficiently.
2. If there is a stain on the surface of the product, remove it immediately using a wet tissue or a wet cloth
3. Store the product neatly after use.
4. Avoid folding or pressing.
5. Damage should be avoided during storage or use.
6. Room temperature storage

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11. Indication
권장사용

Use to protect the body from X-rays as much as possible when working in an environment at risk of radiation exposure.

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Areas that are required of X-ray monitoring should not be covered by this device, otherwise X-ray monitoring could not be implemented.

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13. Labeling
라벨

Inside the bib, you will find the identifies the model name, size, lead equivalency, date of manufacture, site of manufacture, UDI-DI, information of manufacturer

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