INSTRUCTION FOR USE
사용설명서

1. Manufacturer Name
제조사

DR-VU (Co.,Ltd)
(주)디알뷰

Registered no. 5882

2. Manufacturer Address Manufacturing site
주소

40, Dongil-ro 96-gil, Jungnang-gu,Seoul, Republic of Korea
서울시 중랑구 동일로96길40, B1

3. Product name
품목명

Patient positioning or immobilization device
환자고정보조대

4. Product License
품목허가번호

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5. Specification
제품사양

No	Model name	Size(cm)		
		목둘레	폭(앞)	뒤폭(앞기장)
1	Rad-Ban HBA	160	45	2
2	Rad-Ban HBC	160	29	2

6. Intended purpose
사용목적

Medical instrument or machine used to maintain or fix a patient's posture during patient care or injection.

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A. Preparations before use

1. Check the condition of the product before use to check for wrinkles, tearing, and damage.

B. How to use it

1. Place the body-part mat on the treatment table and install the head mat if necessary.
2. Fold the head-part mat in half and place it on the head, wrap the head according to the appropriate length, and attach the Velcro on the head and the Velcro on the body.
3. Attach the shoulder, body, wrist, knee, and leg fixing belt to the Velcro of the body part in the appropriate position, and attach the patient's headband using the Velcro of the back of the head as well.
4. Fixing the shoulder fixing belt in an X shape can increase the fixing power.
5. Attaching the shoulder and body fixing bands so that they overlap slightly can further reduce the patient's movement during treatment.
6. If the patient's legs are thin and the fixing is weak, the fixing force can be strengthened by inserting it between the legs using a doll or cushion and installing Velcro.

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1. It shall not be used for any purpose other than the purpose of use.
2. Be careful not to scratch the surface.
3. Control pollution.
4. Control loss of accessories.
5. When washing, match the Velcro in the same pair so that it does not stick to different parts, and then put it in a laundry net to wash it.

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1. Do not break.
2. In case of contamination, foreign substances such as washing are managed.
3. Manage the loss of the attached fixing belt (by parts).
4. Room temperature storage

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11. Labeling
라벨

Inside the zipper, you will find the identifies the model name, size, lead, date of manufacture, site of manufacture, UDI-DI, information of manufacturer

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